

© Copyright 2001-2008 Synergix Ltd. All rights reserved.

Molecular Conceptor 2

Installation Guide

Standalone version

June 2008

Copyright

© Synergix Ltd. 2001-2008. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language in any form by any means without the written permission of Synergix Ltd.

Software License Notice

Your license agreement that is included with the product specifies the permitted and prohibited uses of the product. Any unauthorized duplication or use of the software in whole or in part, in print, or in any other storage and retrieval system is forbidden.

Licenses and Trademarks

Molecular Conceptor 2 is a trademark of Synergix Ltd.

Windows 98, Windows 2000, Windows ME, Windows XP, Windows Server 2003, Windows Vista, Windows Server 2008 are registered trademarks of Microsoft Corporation.

Novell is a registered trademark of Novell, Inc.

Viewer ActiveX Control is a registered trademark of Accelrys, Inc

HASP is registered trademarks of Aladdin Knowledge System Ltd.

All other trademarks, brands, and product names used in this guide are trademarks of their respective owners.

A. Table of Contents

A.	Table of Contents	3
B.	Molecular Conceptor 2 TM License	4
C.	Introduction	9
C.1	System Requirements	9
C.2	Contents of the Standalone Package	9
C.3	Viewer ActiveX Control	9
C.4	The HASP Key.....	10
D.	Software Installation	11
D.1	Installation of Molecular Conceptor 2	11
E.	Running Molecular Conceptor 2	15
F.	Uninstall procedure	16
F.1	Uninstall Molecular Conceptor 2	16
F.2	Uninstall HASP key Driver.....	16
G.	Questions or Problems	16

B. Molecular Conceptor 2™ License

Synergix®

Synergix Electronic End-User Software License Agreement

THIS AGREEMENT IS A LEGAL DOCUMENT. READ IT CAREFULLY BEFORE COMPLETING THE INSTALLATION PROCESS AND USING THE SOFTWARE. IT PROVIDES A LICENSE TO USE THE SOFTWARE AND CONTAINS WARRANTY INFORMATION AND LIABILITY DISCLAIMERS. BY INSTALLING AND USING THE SOFTWARE, YOU ARE CONFIRMING ACCEPTANCE OF THE SOFTWARE AND AGREEING TO BECOME BOUND BY THE TERMS OF THIS AGREEMENT. IF YOU DO NOT WISH TO DO SO, DO NOT COMPLETE THE INSTALLATION PROCESS. PROMPTLY RETURN THE ENTIRE PACKAGE TO THE PLACE WHERE YOU OBTAINED IT, FOR A FULL REFUND.

1. Definitions

- (a) "Synergix" means Synergix Ltd. and its licensors, if any.
- (b) "Software" means the programs which may also include documentation, associated media (e.g.:text(s), Molecule(s), 3D visualization(s), picture(s), animation(s)), printed materials, voice(s), and online and electronic documentation, supplied by Synergix herewith, and any and all updates thereto.
- (c) "Network License" means a multi users license of the Software accessible across a network controlled by Aladdin License Manager. License terms allow for example 10, 20, 50, 100 (or more) users to access the software concurrently.

2. License

This Agreement allows you to:

- (a) Install and use the Software on a single computer and make one copy of the Software in machine-readable form solely for backup purposes. You must reproduce on any such copy all copyright notices and any other proprietary legends on the original copy of the Software.

(b) If you have purchased a Network License for the Software product, you may install the Software on a storage device, such as a network server, and run the Software on an internal network only, provided the number of concurrent users does not exceed the number authorized by the Network License.

(c) If you have the evaluation version or the demo version of the Software, you may use it for internal testing and evaluation purposes only.

3. License Restrictions

Other than as set in Section 2, you may not install the Software on a storage device, such as a network server, and run the Software on an internal network, you may not use the Software without the supplied protection key, you may not make or distribute copies of the Software, or electronically transfer the Software from one computer to another or over a network. You may not decompile, reverse engineer, disassemble, or otherwise reduce the Software to a human-perceivable form. You may not lend, rent, lease, sublicense, or transfer for value the Software. You may not modify the Software or create derivative works based upon the Software. You may not copy or post any part of it on any network computer or broadcast it in any media. You may not give access to the software from a computer outside your institution, (e.g. access from a web page). You may not transfer the software or any part of it, into any other application, media, or viewer (e.g. PowerPoint, html page). You may not create hard-copy of the software or any part of it.

4. Ownership

The foregoing license gives you limited rights to use the Software. Although you own the media on which the Software is recorded, you do not become the owner of, and Synergix Ltd. retains title to, the Software, and all copies thereof. All rights not specifically granted in this Agreement, including Federal and International Copyrights, are reserved by Synergix.

5. Limited Warranties

(a) Synergix warrants that, for a period of ninety (90) days from the date of delivery (as evidenced by a copy of your receipt): (i) when used with a recommended hardware configuration, the Software and the protection key will perform in substantial conformance with the documentation supplied with the Software; and (ii) that the media on which the Software is furnished will be free from defects in materials and workmanship under normal use. EXCEPT

AS SET FORTH IN THE FOREGOING LIMITED WARRANTY, Synergix DISCLAIMS ALL OTHER WARRANTIES OR CONDITIONS, EITHER EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OR CONDITIONS OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IF APPLICABLE LAW IMPLIES ANY WARRANTIES WITH RESPECT TO THE SOFTWARE AND THE PROTECTION KEY, ALL SUCH WARRANTIES ARE LIMITED IN DURATION TO NINETY (90) DAYS FROM THE DATE OF DELIVERY.

No oral or written information or advice given by Synergix, its dealers, distributors, agents or employees shall create a warranty or in any way increase the scope of this warranty.

(b) (USA only) SOME STATES DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES, SO THE ABOVE EXCLUSION MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER LEGAL RIGHTS WHICH VARY FROM STATE TO STATE.

6. Exclusive Remedy

Your exclusive remedy under Section 5 is to return the Software and the protection key to the place you acquired it, with a copy of your receipt and a description of the problem. Synergix will use reasonable commercial efforts to supply you with a replacement copy of the Software and/or the protection key that substantially conforms to the documentation, provide a replacement for defective media, or refund to you your purchase price for the Software, at its option. Synergix shall have no responsibility if the Software has been altered in any way, if the media has been damaged by accident, abuse or misapplication, or if the failure arises out of use of the Software with other than a recommended hardware configuration. After 90 days, basic cost may apply for replacement of the software and/or protection key.

7. Limitation of Liability.

IN NO EVENT WILL SYNERGIX OR ITS SUPPLIERS BE LIABLE TO YOU FOR ANY CONSEQUENTIAL, INDIRECT, INCIDENTAL, PUNITIVE OR SPECIAL DAMAGES, INCLUDING ANY LOST PROFITS OR LOST SAVINGS, EVEN IF A REPRESENTATIVE OF SYNERGIX OR ANY SUPPLIER HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, OR FOR ANY CLAIM BY ANY THIRD PARTY. SOME

STATES OR JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU.

8. Termination

Without prejudice to any other rights, Synergix may terminate this LA if you fail to comply with the terms and conditions of this LA. In such event, you must uninstall and destroy all your copies of Molecular Conceptor™.

9. General Disclaimer.

ALL CONTENT OF THE SOFTWARE IS PROVIDED "AS IS" AND WITHOUT WARRANTIES OR REPRESENTATIONS OF ANY KIND, EITHER EXPRESSED OR IMPLIED. TO THE FULLEST EXTENT PERMISSIBLE UNDER APPLICABLE LAW, SYNERGIX DISCLAIMS ALL WARRANTIES AND REPRESENTATIONS, EXPRESSED OR IMPLIED, RELATING TO SUCH CONTENT, INCLUDING BUT NOT LIMITED TO ANY WARRANTIES AND REPRESENTATIONS OF NONINFRINGEMENT OF THIRD PARTY INTELLECTUAL PROPERTY RIGHTS, TITLE, MERCHANTABILITY, AND FITNESS FOR A PARTICULAR PURPOSE. ALTHOUGH SYNERGIX WISHES THE CONTENT OF THE SOFTWARE TO BE ACCURATE, COMPLETE, AND CURRENT, SYNERGIX DOES NOT WARRANT OR REPRESENT THAT ANY CONTENT OF THE SOFTWARE IS ERROR-FREE, COMPLETE, OR CURRENT, THAT USE OF SUCH CONTENT WILL BE UNINTERRUPTED, THAT DEFECTS WILL BE CORRECTED. APPLICABLE LAW MAY NOT ALLOW THE EXCLUSION OF CERTAIN WARRANTIES, SO THESE EXCLUSIONS MAY NOT APPLY TO YOU IN WHOLE OR IN PART.

10. (USA only) Government End Users

RESTRICTED RIGHTS LEGEND This software is "Restricted Computer Software." Use, duplication, or disclosure by the U.S. Government is subject to restrictions as set forth in this Agreement and as provided in DFARS 227.7202-1(a) and 227.7202-3(a) (1995), DFARS 252.227-7013 (OCT 1988), FAR 12.212(a)(1995), FAR 52.227-19, or FAR 52.227-14, as applicable.

Manufacturer: Synergix Ltd., Technology Park Malha, Building 1 Jerusalem ISRAEL

11. (Outside of the USA only) Consumer End Users

The limitations or exclusions of warranties and liability contained in this agreement do not affect or prejudice the statutory rights of a consumer; i.e., a person acquiring goods otherwise than in the course of a business.

12. General

This Agreement shall be governed by the internal laws of Israel. This Agreement contains the complete agreement between the parties with respect to the subject matter hereof, and supersedes all prior or contemporaneous agreements or understandings, whether oral or written. All questions concerning this Agreement shall be directed to: Synergix Ltd, Technology Park Malha, Building 1, 91487 Jerusalem ISRAEL Attention: Legal Department.

C. Introduction

Thank you for purchasing *Molecular Conceptor 2*, the unique molecular modeling and drug design courseware.

C.1 System Requirements

The minimum system requirements to install and run *Molecular Conceptor 2* are as follows:

- 350 MHz Intel Pentium processor
- Windows 98SE / Me / 2000 / XP / 2003 / Vista
- 90 MB of free available system RAM
- 720 MB of available disk space
- 24-bit color monitor
- DVD-ROM drive
- Internet Explorer 5.5 or higher

C.2 Contents of the Standalone Package

- HASP key (Hardware protection device)
- *Molecular Conceptor* DVD-ROM
- Client Serial Number Card
- *Molecular Conceptor* Table of Contents
- Installation Guide
- Quick Installation Card

C.3 Viewer ActiveX Control

The 3D visualization of molecules in *Molecular Conceptor 2*™ is based on "Viewer ActiveX Control" from Accelrys Inc

For more information about the Accelrys Viewer products visit <http://www.accelrys.com>

C.4 The HASP Key

The HASP key is a small hardware device (sometimes called a dongle) that protects *Molecular Conceptor 2* against piracy. You must insert the key into your computer only AFTER the completion of the installation procedure.

Molecular Conceptor 2 is distributed either with USB key or Parallel key (the USB key is to be inserted into the USB port, whereas the Parallel key needs to be inserted into the parallel/printer port).



Parallel Key



USB Key

D. Software Installation

D.1 Installation of Molecular Conceptor 2

On Windows 2000/XP/Vista, *Molecular Conceptor 2* must be installed by a user with administrator privileges.

To Install:


1. Close all Windows applications.
2. Insert the *Molecular Conceptor 2* DVD in your DVD-ROM drive. If the "Molecular Conceptor 2 Installation Setup" doesn't automatically run, double-click the 'mc2_setup.exe' program in the root directory of the DVD.




3. Click **Next**.



4. Read carefully the License Agreement. If you accept all the terms, click **I accept the agreement** and press **Next**.




6. Click **Next**.



Regular Installation – Standalone Version: This is the recommended installation for Molecular Conceptor 2 standalone version .

Custom Installation (network version): Not used in Molecular Conceptor standalone version


8. Click **Next**.



The Installation program will copy *Molecular Conceptor 2* files to the computer hard disk, and will install the HASP driver. You might be asked to reboot the computer.

10. Click **Finish**.

11. After the completion of the installation program, insert your HASP key into your computer port (either USB or Parallel)




E. Running Molecular Conceptor 2

Run *Molecular Conceptor 2* via the Start menu:

Start Menu | Programs | Molecular Conceptor 2 | Molecular Conceptor 2

You can also run *Molecular Conceptor 2* by double-clicking on the shortcut icon on your desktop. If you prefer not having the shortcut icon on your Desktop, you can delete it; the program will still be available by the Start menu.



F. Uninstall procedure

During the installation of *Molecular Conceptor 2* two different programs have been installed:

1. Molecular Conceptor 2
3. Hasp Key driver

F.1 Uninstall Molecular Conceptor 2

To uninstall:

1. Open **Control Panel** (e.g. Click **Start Menu** | **Settings** | **Control Panel**).
2. Double-click the **Add/Remove Programs** icon.
3. In the **Add/Remove Programs** dialog box, select **Molecular Conceptor 2**.
4. Click **Remove**.

F.2 Uninstall HASP key Driver

To uninstall:

1. Run the utility program, from the *Molecular Conceptor 2* CD:
'[CD_DRIVE]:\utility\HASP\uninstall_hasp.bat'

G. Questions or Problems

If you have questions or problems with *Molecular Conceptor*, first consult our online help

<http://www.molecular-conceptor.com/support.html>

For additional technical support, contact us at:

support@molecular-conceptor.com